

uMED 20 Vet

Veterinary Defibrillator/Monitor




Saving lives through simplicity



LinkedIn Facebook Instagram

Search Mindray Animal Medical

www.mindrayanimal.com
service@mindrayanimal.com



Fast Defibrillation

Equipped with QShock technology, uMED 20 Vet provide an ultra-fast defibrillation experience. It only takes < 7.5s from power on, charging to ECG recovery. helping to increase the chance of a successful defibrillation.

Intuitive Panel

- Manual defibrillation with clear 1-2-3 steps
- Color-coded buttons for easy identification



Single Operator Design

- External paddles come with buttons for energy selection, charge and shock
- Integrated external paddles supports both large size and small size patients





Higher Efficiency

Advanced Technology

Biphasic defibrillation: improved technology with lower risk of skin burns, less myocardial injury

Impedance Compensation: auto adjustment of peak current and discharge cycles according to the impedance change.

Animal Contact Indication: colored indicators for animal contact condition are available on both main screen and external paddles, help to optimize the contact position or choose the best energy for a better recovery.



Always Standby

Comprehensive Auto-test

Power-off test:

- Daily and weekly test are provided while power off, including the energy delivery test.

Power-on test:

- A comprehensive test whenever the device is turned on, status indicator illuminated once auto test failure.



Status indicator

Durable Battery

- Auto calibration of battery capacity, maintenance free
- Large capacity battery supports longer operation time and life span
- Dedicated for defibrillation application



≥300 shocks @200J